DENTSPLY

ECLIPSE®

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570 West College Avenue
P.O. Box 872
York, PA 17405-0872

www.trubyte.dentsply.com or www.dentsply.com

© 2003, 2004 DENTSPLY International, Inc. All rights reserved.

Printed in U.S.A. 905437 Rev. C 6/08

CE

DeguDent GmbH
Postfach 1364 • 63403 Hanau
GERMANY
Telefon +49/6181/5950

R_xonly

DENTSPLY

ECLIPSE®

PROSTHETIC RESIN SYSTEM

Eclipse Resilient Resin Directions For Use

Indications for Use: Occlusal splints and nightguards.

Contraindications for Use: None known.

Warnings: This product contains polymerizable materials. It does not contain methyl methacrylate (MMA) monomer. It is possible that some dental health care workers may develop an allergy or sensitivity to the product. Wash thoroughly with soap and water after contact. If dermatitis or other symptoms persist, seek medical advice. Do not immerse Eclipse appliances in boiling water.

Severe undercuts need to be blocked out prior to fabrication.

Precautions:

1. This product is intended to be used only as specifically outlined in the Directions For Use. Any use of this product inconsistent with the Directions For Use is at the discretion and sole responsibility of the practitioner.
2. Users with special skin problems, cuts or abrasions may choose to wear protective gloves.
3. When grinding prosthodontic resins, proper ventilation, masks and vacuum systems should be used.
4. Since the Eclipse resins will become the final prosthesis, cleanliness is very important to an acceptable outcome.

Wash hands before using uncured Eclipse materials.

5. Store at 60-80° F (16-27° C), away from moisture and in a light-safe container. The materials may polymerize prematurely if stored at excessively high temperature or if exposed to a visible or UV light source. Do not expose to direct sunlight.
6. Unused Eclipse Resilient Resin should be cured prior to disposal.
7. Items removed from the curing unit are hot! Use heat-resistant gloves provided.
8. Do not remove prostheses from models until cool.
9. If the resin materials will not release easily from the package material (may occur in warm ambient conditions), place the package in a refrigerator/freezer for a few seconds to facilitate removal.
10. The following disinfectants have been used successfully with cured and uncured Eclipse resins:
 - Glutaraldehyde-based disinfectants
 - Quaternary ammonium chlorides
 - Quaternary ammonium chlorides/isopropyl alcohol (low concentration)
 - Spray-based phenolics

Use of other disinfectant products may not be suitable with Eclipse. Avoid disinfectants containing iodine.

Adverse Reactions:

1. Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals.
2. Dust will be generated when grinding these materials. Eye, skin and respiratory irritation may occur if appropriate engineering controls are not used.

NOTE: Stop and read operating manuals for equipment before continuing to read Directions for Use.

DIRECTIONS FOR USE:

1. Paint the model with Al-Cote® Separating Agent and allow it to dry. Separator MUST be dry before fabricating.
2. Place the model in the DENTSPLY® Conditioning Oven to warm, OR warm model with the DENTSPLY® Hot Air Gun for about 1 minute. Model must be warm for Eclipse Resilient Resin to adhere to model. (If you choose to warm the model in the DENTSPLY Conditioning Oven, it is not necessary to bring the model to a temperature of 55°C.)
3. Chill the Eclipse Resilient Resin by running under cold water for 15-30 seconds to facilitate removal from the package.
4. On the warmed model, adapt the Eclipse Resilient Resin based on the design. Severe undercuts need to be blocked out on lingual, buccal and interproximal and a duplicate model made when necessary. Eclipse Resilient Resin usually does not need to be placed below the height of contour of the teeth as with other soft materials. If the model has cooled down, warm it with the DENTSPLY Hot Air Gun and continue to adapt the Eclipse Resilient Resin onto the model. Any excess material can be removed using a sharp instrument.
5. Remove Clear Eclipse Baseplate Resin from package and cut in half length-wise. Adapt the Clear Eclipse Baseplate Resin

directly over the Eclipse Resilient Resin on the occlusal surface. If material and model have cooled down, warm in the DENTSPLY Conditioning Oven or with the DENTSPLY Hot Air Gun to facilitate easier adaptation.

6. Before placing back on the articulator to check occlusion, coat the opposing dentition with Eclipse Model Release Agent. Then, place the model on the articulator and gently tap into occlusion until the pin is down. If needed, additional material can be added. Using the incisal guide pin and table, the nightguard can be put through side-to-side excursions if desired.
7. Smooth the junction of the Clear Eclipse Baseplate and the Resilient Resins with the DENTSPLY® Wax Pencil Pro. Further smoothing can be accomplished with the DENTSPLY Hot Air Gun as well.
8. Remove the model from the articulator and paint the Eclipse material with Eclipse Air Barrier Coating. Immediately process in the Eclipse Processing Unit on Menu #7.
9. Upon completion of the curing cycle: Allow to bench cool to room temperature before rinsing off the Eclipse Air Barrier Coating. Remount on the articulator and verify/adjust occlusion as necessary.
10. Soak splint and model in warm water for 5 minutes to facilitate removal from the model.
11. Wash off Eclipse Air Barrier Coating with water and a soft brush.
12. Trim and smooth borders.
13. Finish and polish.

Optional Method

After Step 4:

- A. Do NOT apply Eclipse Air Barrier Coating. Process the Eclipse Resilient Resin in the Eclipse Processing Unit on Menu #6.
- B. Remove from Eclipse Processing Unit and continue with step 5.

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